covalon
TECHNOLOGIES LTD.

Mobtakeran Salamat Fakher


Exclusive Representative of **Covalon LTD.** in Iran

ColActive[®] PLUS Ag

Collagen Matrix Dressing with Silver

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
ColActive® PLUS

ColActive® PLUS
Collagen Matrix Dressing

WHAT IS IT?

ColActive® Plus is an advanced wound dressing made from collagen, sodium alginate, carboxymethylcellulose (CMC), and ethylenediaminetetraacetic acid (EDTA). **ColActive® Plus** is a topically applied wound dressing that will transform into a soft gel sheet when in contact with wound exudates.

HOW DOES IT WORK?

The dressing maintains a moist wound environment that aids in the formation of granulation tissue and epithelialization. The **EDTA** in the dressing removes zinc to inhibit the activity of matrix metalloproteinases (**MMP**), thus creating a suitable environment for wound healing. **ColActive® Plus** dressings may be trimmed and layered for the management of deep wounds.

WHEN DO I USE IT?

ColActive® Plus is indicated for management of full and partial thickness wounds including: pressure ulcers, diabetic ulcers, ulcers caused by mixed vascular aetiologies, venous ulcers, donor and graft sites, abrasions, traumatic wounds healing by secondary intention, dehisced surgical wounds, first and secondary degree burns.

AVAILABILITY:

COLACTIVE® PLUS

Collagen Matrix Dressing

Reference Number	HCPCS	Size (cm)	Dressings per Carton / per Case
TWBC1016 5550220	A6021 —	5cmx5cm	10/10
TWBC1018 5550440	A6022 —	10cmx10cm	10/10
TWBC1035 5550770	A6023 —	18cmx18cm	10/10

SUGGESTED WOUND TISSUE TYPE:

NECROTIC (After debridement)	✓
SLOUGHY	✓
GRANULATING	✓
EPITHELIALISING	✓
INFECTED	--

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3.

4.

The dressings may be trimmed and layered for the management of deep wounds.

WHEN DO I USE IT?

ColActive® Plus Ag is indicated for management of full and partial thickness wounds including: pressure ulcers, diabetic ulcers, ulcers caused by mixed vascular aetiologies, venous ulcers, donor and graft sites, abrasions, traumatic wounds healing by secondary intention, dehisced surgical wounds, first and secondary degree burns.


AVAILABILITY:

COLACTIVE® PLUS Ag Collagen Matrix Dressing with Silver

Reference Number	HCPCS	Size (cm)	Dressings per Carton / per Case
TWBC1033 5551220	A6021 —	5cm x 5cm	10/10
TWBC1034 5551440	A6022 —	10cm x 10cm	10/10
TWBC1032 5551770	A6023 —	18cm x 18cm	10/10

SUGGESTED WOUND TISSUE TYPE:

NECROTIC (After debridement)	✓
SLOUGHY	✓
GRANULATING	✓
EPITHELIALISING	✓
INFECTED	✓

ColActive® PLUS Ag

Collagen Matrix Dressing with Silver

WHAT IS IT?

ColActive® Plus Ag Collagen Matrix Dressing with **Silver** is an advanced wound care dressing made of collagen, sodium alginate, carboxymethylcellulose (**CMC**), ethylenediaminetetraacetic acid (**EDTA**) and silver chloride.

HOW DOES IT WORK?

The topically applied wound dressing will transform into a soft gel sheet when in contact with wound exudates. The dressing maintains a moist wound environment at the wound surface that aids in the formation of granulation tissue and epithelialization. The **EDTA** in the dressing removes zinc to inhibit the activity of matrix metalloproteinases (**MMP**), thus creating a suitable environment for wound healing. The silver chloride in the dressing provides effective broad spectrum antimicrobial activity.

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Wound Contact Layer

WHAT IS IT?

ColActive® Transfer is pliable, absorbent, and biocompatible, and is designed to intimately contact the wound bed and protect it from secondary dressings (i.e. synthetic foam). ColActive® Transfer has low adherence, which helps minimize wound trauma at dressing changes.

HOW DOES IT WORK?

ColActive® Transfer contains multiple holes that are evenly distributed across the surface area to allow fluid to pass through without impairment, thus maintaining a balanced moist wound environment which has been shown to be beneficial for wound healing. ColActive® Transfer maintains its structural integrity, in various types of wounds, including under negative pressure for a period of up to 5 days (in vitro data).

WHEN DO I USE IT?

ColActive® Transfer is indicated for management of full and partial thickness wounds including pressure ulcers, diabetic ulcers, ulcers caused by mixed vascular etiologies, venous ulcers, second degree burns, donor and graft sites, abrasions, dehisced surgical wounds, and traumatic wounds healing by secondary intention.

AVAILABILITY:



COLACTIVE® TRANSFER

Wound Contact Layer

Reference Number	HCPCS	Size (cm)	Dressings per Carton / Cartons per Case
TWBT1040 5550221	A6021 —	5cm x 5cm	10/10
TWBT1041 5550441	A6022	10cm x 10cm	10/10

SUGGESTED WOUND TISSUE TYPE:

NECROTIC (After debridement)	✓
SLOUGHY	✓
GRANULATING	✓
EPITHELIALISING	✓
INFECTED	--

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Education



- . Kerman University of Medical Sciences, Kerman, Iran - General Practitioner, 2002
- . Mashhad University of Medical Sciences, Mashhad, Iran - Infectious Disease, 2009

Dissertations

- . Determining the relation of diabetes with dysfunction of brain frontal lobe
- . Epidemiological study and some features of the injection drug users hospitalized in the infectious ward of Imam Reza Hospital of Mashhad

Achievements

- . Established Diabetes and Diabetic Foot Clinic in Kerman, Iran - 2018

Dr. Mehdi Khaleghinia
Infectious Disease Specialist

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Courses Attended

- . 3rd Tehran Congress of Hepatitis, Tehran, Iran - 2002
- . 5th Paris Congress of Hepatitis, Paris, France - 2012
- . Good Medical Practice (MEDA), Istanbul, Turkey - 2012
- . Asia Pacific Congress of Hepatitis (APASL), Singapore - 2013
- . International European Congress of Liver and Hepatitis, Amsterdam, Netherlands - 2013
- . 24th Asia Pacific Congress of Hepatitis and Liver, Istanbul, Turkey - 2015
- . 9th Pisa International Diabetic Foot Course (Theory and Practice, Pisa, Italy - 2017
- . 28th International European Congress of Infectious Diseases (ECC MID), Madrid, Spain - 2018
- . International European Congress of Ulcer (EWMA), Krakow, Poland - 2018

Certificates


- . Educational and practical certificate of working with high and low power foot laser for wound
- . Educational and practical certificate of working with ozone therapy for wound

Publications

- . Investigation of the serum levels of IL8, IL6, TGFβETA, and TNFα, by: Reza Bahram Abadi, Mahmoud Sheikh Fatollahi, Seyed Mohammad Hashemi, Hassan Yousefi, Mehdi Khaleghi Nia, and Reza Bidaki, Oxford Academic, 2017

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Diabetic Foot Ulcer

Patient is a 54 year old male with notable past medical history of uncontrolled diabetes. **25** days prior to presentation the patient experienced extensive ulceration involving the phalanxes, metatarsals and forefoot on the right side. Upon presentation **25** days later he was found to have extensive necrosis and infection in the aforementioned sites. Ankle Brachial Index (**ABI**) measurements of the affected limb were noted to be within normal limits. The initial approach consisted of aggressive debridement until viable margins were visualized. Thereafter **ColActive Plus Ag** dressing was applied. The dressing was changed twice per week for the first two weeks for optimal control of infection risk and wound discharge. The dressing was changed weekly thereafter. Treatment was started on **2019-12-01** and concluded on **2020-01-08** with complete epithelization and closure of the lesions.

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
Diabetic Foot Ulcer

A 59 year old female patient with diabetes and a most recent HBA1c measurement of 8.5% presented with a heel ulcer with skin and soft tissue necrosis of 8 months duration. The lesion was highly exudative, producing foul smelling purulent discharge. Examination of the bilateral lower extremities showed a temperature difference in the affected limb and ultrasound studies determined a reduction in blood flow of about 50% to the affected limb. A probe to bone test and scan were negative for osteomyelitis. The patient had a fever of 38 degrees and tenderness over the heel. A diagnosis of necrotizing fasciitis was made and treatment consisted of debridement to uncover viable tissue, followed by application of the collagen matrix contact layer with silver to the wound bed covered with a secondary dressing. Dressing changes were performed three times per week for the first week and the weekly thereafter. Treatment was initiated on 2019-10-11 and concluded on 2019-11-23 with complete healing.

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
2019-01-06

Diabetic Foot Ulcer

A 68 year old female patient with a history of diabetes and undergoing dialysis for end stage chronic kidney disease presented with necrotic fasciitis. The patient underwent surgical fasciotomy and extensive debridement prior to initiating treatment with the collagen matrix contact later with silver. Dressing changes were initially performed twice per week, then weekly thereafter. Treatment was initiated on 2018-12-19 and complete closure was achieved to resolve the case on 2019-01-06.

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
Trauma

A 28 year old male presented with an infected traumatic wound on the right elbow resulting from an injury which occurred approximately 3 months prior. Initial examination revealed necrotic tissues and swab cultures revealed *Pseudomonas aeruginosa*. Following debridement, the wound was dressed with the collagen matrix contact layer with silver and covered with a secondary dressing. For the first 3 weeks of treatment, the contact layer was changed twice per week, followed by weekly changes thereafter. Treatment was initiated on 2019-10-30 and concluded on 2019-12-19 with complete epithelialization and wound closure.

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

Diabetic Ulcer

A 68 year old diabetic male developed and infection of the tissue overlying the tibial region of his forearm where a graft was recently harvested. In addition to history of diabetes, and a last HgA1c of 8%, the patient previously underwent coronary artery bypass grafting. Given the extensive risk factors for vasculopathy, ultrasound doppler performed and identified arterial insufficiency in the affected limb. Following debridement, the wound was dressed with the collagen matrix wound contact layer containing silver. The dressing was changed twice weekly for the first three weeks and then every week thereafter. Treatment was initiated on 2019-10-17 and successfully concluded on 2020-01-02 with complete resolution and would closure.

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



Diabetic Foot Ulcer

A 63 year old male with history of uncontrolled diabetes presented with a non-healing ulcer on the left foot following mid-metatarsal amputation. No granulation was evident at was evident 25 days after the surgery, at which point the patient started to develop infection of the site making below the knee amputation a possibility. Following aggressive antibiotic treatment and debridement of the infected wound, topical treatment with a hydrogel and ozone therapy was initiated. The collagen matrix contact layer with silver was applied and initially changed twice per week for 4 weeks then once weekly thereafter. Treatment was initiated on 2019-6-21 and concluded on 2019-08-21 with successful closure of the lesion.

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

2020-01-08

Diabetic Foot Ulcer

A 58 year old female presented with ulcer in the calcaneus region of the right heel. The patient had a history of uncontrolled diabetes (last HgA1c 9.5%) and arterial insufficiency. As anticipated, the doppler ultrasound imaging demonstrated narrowing in the peroneal region. She did not consent to stenting of the lesion, so it was not pursued despite the patient being a good candidate for the procedure. Autolytic debridement of the lesion was carried out and the wound was dressed with the collagen matrix wound contact later with silver. The dressing was initially changed twice per week for three weeks and once per week after that. Treatment was started on 2019-10-16 until complete resolution on 2020-01-08 .

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
Trauma

A 5 year old boy initially presented with a recent and extensive trauma to the right hand. He was subsequently lost to follow up and later presented with necrosis of the wound. After sharp debridement followed by hydrogel and ozone therapy, the wound was dressed with the collagen wound contact later with silver. The dressing was changed twice per week for the first 4 weeks and weekly thereafter. Treatment was started on 2018-09-21 and concluded on 2018-10-18 with complete resolution.

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
Diabetic Foot Ulcer

A 63 year old female with necrotizing fasciitis on the plantar surface of the right foot involving the toes. The patient had extensive medical history including hepatitis B, Liver cirrhosis and uncontrolled diabetes (HbA1c 10.5%). Following intensive antibiotic therapy, the involved were amputated and the collagen matric contact layer with silver was applied. Dressing changes were performed three times per week for the first two weeks followed by once per week. Treatment was started on 2018-04-23 and concluded on 2018-11-11 with epithelization and closure of the lesion.

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

2019-11-08

Diabetic Foot Ulcer

A 63 year old female with a history of diabetes (HgA1c of 9%) initially presented with the chief complaint of pain in the region of the right foot. The patient believed in traditional healing methods such as using leeches. Examination revealed three lesions on her lower extremities and necrotic fasciitis of the 4th and 5th metatarsals, as well as a lesion involving the upper calf region. She was also found to have arterial insufficiency of the iliac artery and was successfully stented. The lesions were debrided and dressed with the collagen wound contact layer containing silver, which were changed twice weekly for the first two weeks and weekly thereafter. Treatment was started on 2019-08-09 and finished on 2019-11-08 with resolution.

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
2020-01-16

Leishmaniasis of the forearm

A 65 year old diabetic patient (HgA1c of 7.5%) presented with Leishmaniasis of the forearm. She underwent medical treatment for two months. She then presented with a raised and infected lesion at the affected site. The site was found to be infected with Staph Aureus. She then underwent photodynamic therapy along with application of the collagen matrix contact later containing silver. Initially, dressing changed occurred twice per week and then weekly thereafter. Application of the collagen matrix contact layer was started on 2019-11-16 until complete resolution of the symptoms on 2020-01-16.

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Sales Office :

Kashani Expy, Noor Building,
No.464, 2nd Fl. Unit 9, Tehran,
Iran

Web & E-mail

www.m-s-f.net
E-mail : info@m-s-f.net

Phone & Fax

+982144086988
+982186081472

Address

Farhangian Ave, Moalem Blv,
Yas Medical Complex, 5th
Fl. RADIN Diabetes &
Diabetic Foot Clinic.
Kerman, Iran

Web

www.drkhaleghinia.com
E-mail : info@drkhaleghinia.com

Phone & Fax

+983432127407

Design by : Amir Farhad Fakhredin

MOBTAKERAN SALAMAT FAKHER Department of Design and Graphics